



ADVANCED NUCLEAR CONSULTANTS



Philips Cardio 60/ Cardio EPIC Gamma Camera

The Philips Cardio 60 and the Cardio EPIC are dedicated cardiac cameras featuring an open, patient friendly design. Both cameras are capable of imaging myocardial perfusion studies in both the prone and supine positions. Both cameras have the capability to perform general nuclear imaging studies as well.

Features & Benefits

- Fixed 90 degree detectors
- Semi-automatic collimator exchange
- Optimized 90° SPECT detector configuration
- Rectangular 200mm x 150mm FOV

ANC's Standard Equipment

- Uninterruptible Power Supply (UPS)
- Network Printer
- Cardiac IVY Gate
- Bar Phantom

ANC's Standard Services

- Equipment Applications Training
- 24 Hours on Site Response Time
- Less Than One Hour Phone Response
- Unlimited Phone Support
- Free Equipment Shipping and Freight Charges

Options to Consider

- Hot Lab and Ancillary Equipment
- Patient Comfort Kit
- Extended Service Agreement
- Preventive Maintenance
- Archive Tracker (Data Back up)

Thank you very much for your inquiry. We appreciate the opportunity to serve you!